



HEALING ZONE AESTHETICS™

a division of Diamond Aesthetic Systems

12 Point Face & Neck Contouring & Rejuvenating Program

Maintaining a toned vital complexion requires a total full face approach to aesthetics. Areas of the face have different thickness of skin, supporting bone structure and muscle configuration. Being specific to the needs of each area of face, recognizing the environmental influences and contributing stress factors, allows for aesthetic programs to deliver optimal results. Add the Neck and Decollette for a complete treatment approach.

***Collajam & Cabroxy Power Face Masks** are a great complement to all skin care therapies.

Equipment: Microdermabrasion, SmoothSculpt, Micro-current, EMS, Oxygen Infusion, RAF

MesoTherapy Product: Dermal Natural Performance, CoEnzyme Q10, DMAE, Growth Factor, Muscle Striz, Brite White, Tri-Tide, Muslce Striz, O24U, C-Plus

Area of Focus	Equipment	MesoTherapy with Large Molecule Electroporation/ Dermal Pen or MicroRoller	Purpose
Deep Lines & Wrinkles	Microdermabrasion/RAF	DMAE	To build collagen, refine surface, soften lines & plump dermis
Heavy eyelids & eyebrows	SmoothSculpt/ Micro-current	DMAE	To break-up congestion and stagnation that weight down skin and then provide eye-lift
Puffy undereye pouches	SmoothSculpt/ Vacuum tech	TriTide/CoEnzyme Q10	to remove fluids and toxins That collect under the skin and cause stagnation and skin stretching
Dark Circles	SmoothSculpt/ Vacuum tech	Brite White & Dermal Natural Performance	lessen and clear stagnation under the eyes that contributes to darkness/clear surface pigmentation
Lack of Cheekbone Definition	SmoothSculpt/ Micro-Current/ RAF	Dermal Natural/ Muscle Striz	remove stagnation and old fat that has settled under the cheekbones and then tone and tighten skin

Area of Focus	Equipment	MesoTherapy with Large Molecule Electroporation/ Dermal Pen or MicroRoller	Purpose
Droopy jowls and lack of jawline Definition	SmoothSculpt/ EMS Technology/ RAF	Muscle Striz/Growth Factor/Dermal Performance	to remove fat accumulation in sagging skin, to tighten and tone skin, to stimulate collagen, restore muscle strength
Pigmentation	Light microderm/O2 Infusion/ topicals	Brite White/O24U	To suppress melanin and lessen dark spots and even skin tone
Acne/tempermental breakouts	SmoothSculpt/Microderm/ O2 Infusion	O24U/DMAE/Dermal Natural	Clear, kill and reduce dermal bacteria, stimulate cleansing lymph, oxygenate skin tissue, freshen and strengthen skin surface
Rosacea	O2 Infusion	DMAE /O24U/ Dermal Natural	Reduce inflammation, provide healthy balanced blood flow to clear skin irritation
Dehydrated skin With fine lines	Microdermabrasion/ O2 Infusion	Growth Factor/Co Enzyme Q-10/ Brite White	to feed and nourish skin cells, refine skin surface to allow for better penetration of product, to tighten and rejuvenate skin
Neck and Under Chin	SmoothSculpt/Microderm/EMS Technology	Growth Factor/DMAE/ Muscle Striz	to reduce wrinkles, smooth tone and tighten the skin, reduce excess fat and fluids under the chin (double chin) strengthen muscles to provide lift and lessen folds at bottom of the neck
Decollette	Microderm/O2 Infusion	Brite White, CoEnzyme Q 10	Smoothing and Refining Skin lines and wrinkles on the Surface and to brighten and lighten the skin



ARM TONING & SKIN NOURISHING PROGRAM

This program is designed to smooth and restore the contour of the arms. Smoothing out bumps and lumps, nourishing the skin surface while building tone and definition to the muscle are all features of this program.

*Carboxy Gel Power Body Masks are a great complement and support to the Arm Toning & Skin Nourishing Program.

Equipment: SmoothSculpt/Vacuum Technology, Electro Muscle Stim, O2 Infusion, Microdermabrasion

MesoTherapy Product: Lipoactivex – Ampeloredux III - Muscle Striz – Biocomplex Oligorevitalizing

Area of Focus	Equipment	MesoTherapy with Large Molecule Electroporation/ Dermal Pen or MicroRoller	Purpose
Upper Arm: Loose Flabby skin	SmoothSculpt/Vacuum Technology	Lipoactivex/Ampeloredux III	reduce fat & fluid filled tissue that causes droopy skin
Upper Arm: Lack of Muscle Contour	Electro Muscle Stim	Reaffir Si/Muscle Striz	restore toned definition
Back of Arms/rough bumpy skin	Microdermabrasion/ O2 Infusion	Biocomplex Oligorevitalizing	restore smooth healthy skin
Elbows/dry skin folds	Microdermabrasion	Biocomplex Oligorevitalizing	smooth and revitalize skin
Forearm/dehydrated thin skin	Microdermabrasion/ O2 Infusion	Growth Factor/Biocomplex	build collagen/nourish skin
Hands (brown spots)	Microdermabrasion	Growth Factor/Co Enzyme Q-10/ Brite White	build collagen/nourish skin Fade and retard melanin



ABDOMEN DISTENSION & FAT REDUCTION PROGRAM

Whether this condition is caused by excess belly fat, lack of muscle tone, poor organ strength, improper digestion or stress the abdomen is the center focus of the body. It needs to be strong, healthy and toned. Clothes fit better, posture improves and a general sense of well being prevails when our core muscles are working.

*Carboxy Gel Power Body Masks are an excellent complement to all body toning treatments.

Equipment: : SmoothSculpt/Vacuum Technology- Eletro Muscle Stimulation – Radio Frequency

MesoTherapy Product: Lipoactivex – Ampeloredux III - Muscle Striz – Biocomplex Oligorevitalizing

Area of Focus	Equipment	MesoTherapy with Large Molecule Electroporation/ Dermal Pen or MicroRoller	Purpose
Waist Roll	SmoothSculpt/Vacuum Technology Radio Frequency	Lipoactivex/Ampeloredux III	To move stagnation & breakdown stored fat while tightening muscles and provide contours to the figure
Hip Fat Pockets	SmoothSculpt/Vacuum Technology	Reaffir Si/Muscle Striz	To release stored fats and fluids
Belly Fat	SmoothSculpt/Vacuum Technology Electro Muscle Stimulation Radio Frequency	Ampeloredux III/ Muscle Striz/ Biocomplex Oligorevitalizing	Moving & breaking down old stored fat. Tightening muscles & restoring firmness. Improving skin tone and elasticity. Stimulate better digestion/organ health.



HEALING ZONE AESTHETICS™

a division of Diamond Aesthetic Systems

LEG CONTOURING & SHAPING PROGRAM

This program is designed to smooth the skin and contour of the leg by moving areas of stagnation, clearing fat pockets, stimulating the lymph system to remove excess fluids and debris as well as to strengthen and tone muscles.

*Carboxy Gel Power Body Masks are an excellent support treatment for the Leg Contouring and Shaping Program

Equipment: SmoothSculpt, ElectroMuscle Stimulation/Electroporation

MesoTherapy Product: Lipoactivex – Ampeloredux III – Reaffir Si - Muscle Striz

Area of Focus	Equipment	MesoTherapy with Large Molecule Electroporation/ Dermal Pen or MicroRoller	Purpose
Outer Thighs/Cellulite “mid thigh to hip”	SmoothSculpt/Vacuum Technology	Lipoactivex/Ampeloredux III	New & Old Fat Breakdown Cleansing Lymph Drainage
Inner Thighs	SmoothSculpt/EMS	Reaffir Si/Muscle Striz	To tighten & tone skin tissue & strengthen muscle
Back of Thighs	EMS/Smooth Sculpt/Vacuum	Reaffir Si/Muscle Striz	To smooth ridges and tighten muscles
Knees & Ankles	Lymph Stimulating Equipment	Lipoactivex	Stimulating lymph drainage for removal of excess fluids
Buttock	Smooth Sculpt	Ampeloredux III	To breakdown stored fats that weigh down the buttock causing droop and ridges to the back of thighs